
TR ARTERIAL KAN GAZI ANALİZİ İÇİN İĞNELİ
ŞİRINGA KULLANIM TALİMATI

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GB SYRINGE FOR ARTERIAL BLOOD GAS ANALYSIS INSTRUCTION FOR USE

Product Type

Synges for Blood Gas Analysis	Hypodermic Needle	
Luer 3 pieces Syringe 2 ml	22G	

Intended Use

Synges for Arterial Blood Gas Analysis are made of empty, sterile and disposable plastic materials, are used by end users for the blood aspiration. The syringes are primarily for human use.

Synges for Arterial Blood Gas Analysis are suitable for use with luer tip connection connectors in accordance with EN ISO 80369-7.

Sterile single use syringes for Arterial Blood Gas Analysis are not intended for injection, but are used to withdraw blood to be used in gas analysis by using a hypodermic needle. The anticoagulant solution, lithium heparin, is sprayed into the syringe. The anticoagulant prevents the blood coagulation.


Product Components

Name	Raw material
Barrel	PP – Polypropylene
Plunger	PP – Polypropylene
Gasket	Isoprene rubber, Latex free.
Hypodermic Needle	Stainless Steel, SS304, PP – Polypropylene
Anticoagulant	80 IU Liquid Heparin Solution
Cap	PP – Polipropilen
Packaging	Medical grade paper and transparent film for EtO (individually blister packed)

Instruction for Use

- Open the package by peeling it from the place marked by the arrow.
- Visually inspect that the content of package is complete and the package is undamaged.
- The product packaging includes the syringe and needle together.
- During priming, hold the syringe with one hand and the needle with the other.
- Make sure that the grading lines on the syringe cylinder are legible.
- Without removing the protective cap of the needle, firmly fix the syringe and needle together by pushing the needle towards the syringe tip and turning it clockwise.
- Carefully remove the needle's protective cap without damaging the needle.
- After sampling, immediately evacuate air bubbles from the syringe on to a gauze swab and close the cap. If liquid is to be injected into the body with a syringe; draw the required amount of liquid into the syringe.

ARTERIAL KAN GAZI ANALİZİ İÇİN İĞNELİ ŞIRINGA KULLANIM TALİMATI / SYRINGE WITH NEEDLE FOR BLOOD GAS ANALYSIS INSTRUCTION FOR USE

Dok./Doc.No: **GKP.46** Tarih/Date: 1.10.20205 Rev.No: 0

Sayfa No/Page No: 7 /9

- After removal of air rotate the syringe around itself to homogenize blood and heparin.
- According to the recommendations from CLSI, the time between drawing blood and performing analysis should not exceed 30 minutes at room temperature.
- Before analysis, the dead space in the syringe must be expelled by flushing a drop of blood into a pad.

Intended Patient Population

There is no limitation for the patient population.

Intended User

Healthcare personnel

Device Description Intended for Use With The Device

Suitable for use with connectors "EN ISO 80369-7 Small diameter connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravenous or hypodermic applications".

Indications


Blood gas analysis provides information about patients' respiratory and metabolic status, with measurements of pH, carbon dioxide (pCO₂) and oxygen (pO₂), as well as electrolytes, lactate and hemoglobin. Syringes for Arterial Blood Gas Analysis are used to collect blood to be used in gas analysis. The anticoagulant it contains prevents blood from clotting. Approximately 80 IU of liquid heparin solution (porcine intestinal mucosa) is applied to the syringe, and the heparin solution is subsequently dried in a drying oven. The product shelf life is 2 years.

Contraindications

Syringes for Blood Gas Analysis should not be used in more than one patient. The product can be used in a one-time sterile condition with a single patient. Cleaning or resterilization is not possible. It is not intended for use other than its stated intended use. Dispose of as medical waste after use.

Cautions

- Do not use the product if the package is damaged.
- Visually inspect that the content of package is complete and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- Before use, check that the hypodermic needle is not blocked and the flow is not obstructed.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused or re-sterilized.
- Do not touch the needle tip to prevent contamination/ needle sticks. Skin puncture with a contaminated needle may result in serious illness such as hepatitis, AIDS and known/unknown diseases.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.
- Work within the specified temperature and time between blood collection and analysis.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.

ARTERIAL KAN GAZI ANALİZİ İÇİN İĞNELİ ŞİRINGA KULLANIM TALİMATI / SYRINGE WITH NEEDLE FOR BLOOD GAS ANALYSIS INSTRUCTION FOR USE

Dok./Doc.No: **GKP.46** Tarih/Date: 1.10.20205 Rev.No: 0

Sayfa No/Page No: 8 /9

Hypodermic Needle Configurations

Needle Configurations	
Gauge	22G
Color	Black
Needle Length	1 1/4"

Storage Conditions

- Store at an ambient temperature of 2-32°C.
- Keep it away from direct sunlight and moist environment.
- Place and keep the box in the direction of arrow.
- Do not place more than 5 successive boxes on the top of each other.
- The shelf life of product is 2 years.

Description of Symbols Used

	Manufacturer		Fragile. Handle with care
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Distributor		Move vertically
	Batch code		Temperature limit 2-32°C
	Catalogue number		Humidity limitation (20% – 80%)
	Medical device		Do not use if package is damaged and consult instruction for use

	Sterilized using ethylene oxide		Do not reuse
	Do not resterilize		Latex free
	Single sterile barrier system		Non-pyrogenic
	Single sterile barrier system with protective packaging outside		Consult instruction for use
	Marking for meeting the requirements of all relevant European Medical Device Directives/European Legislation The number displayed on the bottom is the Notified Body Number		Caution

Manufacturing Site

Set Medikal Sanayi ve Ticaret A.Ş.
Osmangazi Mah. Mareşal Fevzi Çakmak Cad. No:18 Esenyurt – İstanbul Türkiye
Tel: 0212 622 04 00
www.setmedikal.com.tr

Distributor

NT MEDİKAL
Necati TÜRKMEN
YEŞİLOBA MAHALLESİ 46205 SOKAK NO:6/D
Seyhan Adana Türkiye
0533 647 5023
www.ntmedikal.com info@ntmedikal.com